

Calibration of a GC System for the Quantification of CO₂ Capture in Aqueous Solutions of Alcanolamine

Alain Valtz Dr. Christophe Coquelet, CEP (Center of Energy and Processes) Mines ParisTech, France

<http://www.mines-paristech.eu/Researcher/Fields-of-Research/Energy-and-process-engineering/CEP/>

INTRODUCTION

MINES ParisTech is one of France's oldest higher education institutions. While its inspirational model has not changed over the years, its teaching methods and research have been at the leading edge since 1783. Two Nobel laureates trained at the School: Georges Charpak (Physics, 1992) and Maurice Allais (Economics, 1988).

MINES ParisTech is established in five major fields: Earth sciences and the environment; Energy and process engineering; Mechanical engineering and materials; Mathematics and systems; Economics, management and society.

The Center of Energy and Processes (CEP) has skills useful for studying the transformation of matter and energy in many fields, focussing on **complex energy systems**, particularly in transient conditions, and on **controlling their emissions**.

The TEP laboratory specializes in the measurement of thermodynamic properties http://www.fcep.mines-paristech.fr/Accueil_Eng.htm. This laboratory, "Thermodynamique des équilibres entre phases", provides industry with specialists skilled in both experimental and theoretical aspects of handling fundamental and applied works successfully.

DESCRIPTION

A key area of study for the TEP lab is the measurement of partition coefficients (ELV) of mercaptans for aqueous amines in the presence of carbon dioxide and/or hydrogen sulfide.

The release of CO₂ is measured using a GC method with a TCD detector. A calibration regime is essential because of the quantitative method of analysis. Samples may be either liquid or gaseous.

Calibrations are required before, during and after the testing process which takes place over several months. The accuracy and precision of the calibrations are critical.

The GC system is calibrated by injecting a known volume of CO₂ aspirated directly from a gas bottle through a septum. This is done with manual gas tight syringes and is operator dependent, often resulting in inaccurate measurements. To reduce discrepancies a magnifying tool was developed to improve the accuracy, precision and consistency of syringe use by different operators. Even with the new tool the calibration process remains time consuming and operator dependent.



Figure 1. Magnifying tool developed to improve the accuracy, precision and consistency of syringe use.

To address these methodology issues, CEP/TEP lab have evaluated and now use eVol[®], the world's first digital analytical syringe.

METHOD

eVol[®] is used to collect CO₂ from a gas bottle through a regulator equipped with a septum port. Using the 50 µL and 500 µL eVol[®] XCHANGE[®] syringes, volumes of gas are aspirated and injected into the GC system injection port. The CO₂ concentration is directly proportional to the sample volume. It elutes on a packed column in less than 50 seconds with levels detected by a TCD detector.

Each time the total aspirated volume is injected. Increasing volumes of aspirated sample are successively injected to set up the calibration curve with the test for each volume repeated several times.

Tests are performed in a stable and controlled environment at 23 °C and an atmospheric pressure of 1006.3 mBar.

RESULTS

Calibration using eVol® 50 µL syringe

Injection of 64 aspirated volumes from 2 µL to 50 µL.




Figure 2. Linear regression of the calibration curve, $R^2 = 0.9991$




Figure 3. Accuracy of several injected volumes from 2 µL to 50 µL.

Calibration using the 500 µL syringe

Injection of 82 aspirated volumes from 20 µL to 500 µL.




Figure 4. Linear regression of the calibration curve.




Figure 5. Accuracy of several injected volumes from 20 µL to 500 µL.

DISCUSSION

The accuracy of low volumes (2.00 µL and 2.5 µL) using the eVol® 50 µL syringe is affected by Helium back pressure, which causes a dilution of CO₂ during the injection process.

The dispersion observed for low volumes using the 50 µL syringe and reported in figure 2 is not observed in figure 4 when the 500 µL syringe is used.

CONCLUSIONS

Using eVol® for injection of different volumes of CO₂, CEP/TEP was able to achieve calibration of their instrument with levels of accuracy and precision not previously obtained using manual syringes.

In addition, the needles on the eVol® syringes allow for sample collection directly from the gas bottle – similarly to the manual syringes used previously.

Using the eVol®, CEP/TEP lab increased the accuracy of their calibration while reducing the handling time involved in this process.

For more information contact our technical customer support team on: techsupport@sge.com

AUSTRALIA & PACIFIC REGION

SGE Analytical Science Pty Ltd
Toll Free: 1800 800 167
Tel: +61 (0) 3 9837 4200
Fax: +61 (0) 3 9874 5672
Email: support@sge.com

CHINA

SGE Shanghai Representative Office
Tel: +86 21 6407 9382
Fax: +86 21 6407 9386
Email: china@sge.com

MIDDLE EAST

SGE Gulf
Tel: +971 6 557 3341
Fax: +971 6 557 3541
Email: gulfsupport@sge.com

EUROPE

SGE Europe Ltd
European Head Office
Toll Free: 00800 2790 8999
Toll Free Fax: 00800 2626 2609
Tel: +44 1908 568 844
Fax: +44 1908 566 790
Tel France: +33 1 69 29 80 90
Fax France: +33 1 69 29 09 25
Tel Germany: +49 (0) 6155 / 60746 0
Fax Germany: +49 (0) 6155 / 60746 50
Email: europe@sge.com

INDIA

SGE Laboratory Accessories Pvt Ltd
Tel: +91 22 24715896
Fax: +91 22 24716592
Email: sgeindia@vsnl.com

UNITED STATES OF AMERICA

SGE Incorporated
Toll Free: (800) 945 6154
Tel: +1 512 837 7190
Fax: +1 512 836 9159
Email: usa@sge.com

JAPAN

SGE Japan Inc
Tel: +81 45 222 2885
Fax: +81 45 222 2887
Email: japan@sge.com

