

Syringes | Guided Plunger

Diamond Syringe - Syringe Brilliance

Guided plunger syringes are the most rugged syringe available, making them perfect for industrial environments.

Features and Benefits

- Extended barrel guides plunger during injection.
- Robust and rugged.

Recommended Applications

Ideal for student and industrial use.



Product Specifications

Available in 5 μ L to 10 μ L capacities.

Accuracy and Reproducibility	± 1 % (displaced volume)
Temperature Range	5-70 $^{\circ}$ C Fixed Needle, 5-120 $^{\circ}$ C Removable Needle
International Standards Traceability	

For more information about this product visit www.sge.com or contact techsupport@sge.com.