

Injection Port | Signet Septa MN Grade

Give Your GC the Seal of Approval

Injection Port | Septum

Desired septa attributes:

- Reliably seal against the carrier gas pressure in the inlet.
- Capable of being pierced and resealing time after time without pieces of the septa being deposited in the GC inlet system.
- Does not contaminate or bleed material into the chromatographic system.



Features and Benefits

- Premium septum for autosamplers with up to 400 injections per septum.
- Comes supplied in a blister pack which prevent absorption from the environment.

Recommended Applications

Premium septa for autosamplers with up to 400 injections per septum.

Product Specifications

Sizes available (mm): 5, 9, 11, 11.5, 12.5, 17

Material	Durability	Resealing	Solvent Resistance	Tear Resistance	Maximum Injector Temperature
High temperature silicone	Excellent	Excellent	Excellent	Excellent	350 °C, 300 °C for 17 mm size

For more information about this product visit www.sge.com or contact techsupport@sge.com.